

Item Catalogue

NUPCO CODE: 4232150044000 | SN: 1231

NUPCO ITEM DESCRIPTION

STEM, LONG FOR RADIAL HEAD REPLACEMENT, ASSORTED SIZES AS REQUESTED BY END USER

VENDOR PRODUCT DESCRIPTION

No TR-SL long-stem SKU family is listed in the specified Anatomic Radial Head catalogue; source lists sterile stem implants.

SKU CODES INCLUDED FOR THIS NUPCO CODE

- TR-S0600-S
- TR-S0602-S
- TR-S0604-S
- TR-S0608-S
- TR-S0700-S
- TR-S0702-S
- TR-S0704-S
- TR-S0708-S
- TR-S0800-S
- TR-S0802-S
- TR-S0804-S
- TR-S0808-S
- TR-S0900-S
- TR-S0902-S
- TR-S0904-S
- TR-S0908-S
- TR-S1000-S
- TR-S1002-S
- TR-S1004-S
- TR-S1008-S
- TR-S0606-S
- TR-S0706-S
- TR-S0806-S
- TR-S0906-S
- TR-S1006-S




Surgical Technique



System Features [continued]

Instrumentation



Broaches

The innovative broaches in the Anatomic Radial Head System allow the surgeon to create a precise opening in the radial canal for proper insertion of the implant. The broaches enter the radial canal in a straight direction and are less likely to broach the canal at an angle, thus preventing improper implant placement. Spiral flutes on the broach are designed to displace bone during broaching. The implant stem is 0.5 mm larger than the broach diameter. The trial stem diameter is 0.5 mm smaller than the broach diameter to allow for ease of trial insertion and removal.

A mallet should be used to insert the broach. Side pegs are provided for removal with a mallet and also provide a T-handle to ease broach insertion and extraction. Furthermore, the broaches are color-coded for ease of trial implant selection.



Ordering Information

Tray Components

Instrumentation

1	6.0 Stem Gauge Assembly	TR-TGA06
2	Height Gauge	TR-TG02
3	Trial Gauge	TR-TG01
4	Quick Release T-Handle	MS-T1212
5	Stem Removal Tool	TR-MS30
6	Impactor Base	TR-MS03
7	6 mm Collar Reamer	TR-CRA06
8	7 mm Collar Reamer	TR-CRA07
9	8 mm Collar Reamer	TR-CRA08
10	9 mm Collar Reamer	TR-CRA09
11	10 mm Collar Reamer	TR-CRA10

Trial Head Implants

12	20 mm Trial Head, Left	TR-TH20L
13	20 mm Trial Head, Right	TR-TH20R
14	22 mm Trial Head, Left	TR-TH22L
15	22 mm Trial Head, Right	TR-TH22R
16	24 mm Trial Head, Left	TR-TH24L
17	24 mm Trial Head, Right	TR-TH24R
18	26 mm Trial Head, Left	TR-TH26L
19	26 mm Trial Head, Right	TR-TH26R
20	28 mm Trial Head, Left	TR-TH28L
21	28 mm Trial Head, Right	TR-TH28R

Trial Stem Implants

22	6 mm x 0 mm Trial Stem	TR-TS60
23	6 mm x 2 mm Trial Stem	TR-TS62
24	6 mm x 4 mm Trial Stem	TR-TS64
25	6 mm x 8 mm Trial Stem	TR-TS68
26	7 mm x 0 mm Trial Stem	TR-TS70
27	7 mm x 2 mm Trial Stem	TR-TS72
28	7 mm x 4 mm Trial Stem	TR-TS74
29	7 mm x 8 mm Trial Stem	TR-TS78
30	8 mm x 0 mm Trial Stem	TR-TS80
31	8 mm x 2 mm Trial Stem	TR-TS82
32	8 mm x 4 mm Trial Stem	TR-TS84
33	8 mm x 8 mm Trial Stem	TR-TS88
34	9 mm x 0 mm Trial Stem	TR-TS90
35	9 mm x 2 mm Trial Stem	TR-TS92
36	9 mm x 4 mm Trial Stem	TR-TS94
37	9 mm x 8 mm Trial Stem	TR-TS98
38	10 mm x 0 mm Trial Stem	TR-TS100
39	10 mm x 2 mm Trial Stem	TR-TS102
40	10 mm x 4 mm Trial Stem	TR-TS104
41	10 mm x 8 mm Trial Stem	TR-TS108

Ordering Information [continued]

Sterile Components

Head Implants

20.0 mm Head, Right	TR-H200R-S
20.0 mm Head, Left	TR-H200L-S
22.0 mm Head, Right	TR-H220R-S
22.0 mm Head, Left	TR-H220L-S
24.0 mm Head, Right	TR-H240R-S
24.0 mm Head, Left	TR-H240L-S
26.0 mm Head, Right	TR-H260R-S
26.0 mm Head, Left	TR-H260L-S
28.0 mm Head, Right	TR-H280R-S
28.0 mm Head, Left	TR-H280L-S

Stem Implants

6.0 mm x 0.0 mm Stem	TR-S0600-S
6.0 mm x 2.0 mm Stem	TR-S0602-S
6.0 mm x 4.0 mm Stem	TR-S0604-S
6.0 mm x 8.0 mm Stem	TR-S0608-S
7.0 mm x 0.0 mm Stem	TR-S0700-S
7.0 mm x 2.0 mm Stem	TR-S0702-S
7.0 mm x 4.0 mm Stem	TR-S0704-S
7.0 mm x 8.0 mm Stem	TR-S0708-S
8.0 mm x 0.0 mm Stem	TR-S0800-S
8.0 mm x 2.0 mm Stem	TR-S0802-S
8.0 mm x 4.0 mm Stem	TR-S0804-S
8.0 mm x 8.0 mm Stem	TR-S0808-S
9.0 mm x 0.0 mm Stem	TR-S0900-S
9.0 mm x 2.0 mm Stem	TR-S0902-S
9.0 mm x 4.0 mm Stem	TR-S0904-S
9.0 mm x 8.0 mm Stem	TR-S0908-S
10.0 mm x 0.0 mm Stem	TR-S1000-S
10.0 mm x 2.0 mm Stem	TR-S1002-S
10.0 mm x 4.0 mm Stem	TR-S1004-S
10.0 mm x 8.0 mm Stem	TR-S1008-S

Instruments

Stem Removal Tool	TR-MS30
Total Radial Head Tray Lid	TR-0002
Total Radial Head Tray Base	TR-0003
Total Radial Head Top Tray	TR-0004
Total Radial Head Bottom Tray	TR-0005
Total Radial Head Trial Caddy	TR-0006

Optional Items

6.0 mm x 6.0 mm Stem	TR-S0606-S	8.0 mm x 6.0 mm Trial Stem	TR-TS86
7.0 mm x 6.0 mm Stem	TR-S0706-S	9.0 mm x 6.0 mm Trial Stem	TR-TS96
8.0 mm x 6.0 mm Stem	TR-S0806-S	10.0 mm x 6.0 mm Trial Stem	TR-TS106
9.0 mm x 6.0 mm Stem	TR-S0906-S	ARH STD & Optional Trial Gauge	80-0832
10.0 mm x 6.0 mm Stem	TR-S1006-S	ARH STD & Optional Trial Caddy Base	80-0833
6.0 mm x 6.0 mm Trial Stem	TR-TS66	ARH STD & Optional Trial Caddy Lid	80-0857
7.0 mm x 6.0 mm Trial Stem	TR-TS76		