

Item Catalogue

NUPCO CODE: 4232191207900 | SN: 1232

NUPCO ITEM DESCRIPTION

PLATE, VARIABLE ANGLED LCP CLAVICLE PLATES 2.7 LATERAL ASSORTED SIZES, SHAPES, AND TYPES, AS REQUESTED BY END USER

VENDOR PRODUCT DESCRIPTION


Clavicle Plating System superior distal/lateral extension clavicle plates

SKU CODES INCLUDED FOR THIS NUPCO CODE

- 70-0111
- 70-0112
- 70-0116
- 70-0117
- 70-0123
- 70-0124
- 70-0125
- 70-0126
- 70-0319
- 70-0320

Superior Distal Clavicle Plate Surgical Technique [continued]

Figure 8




6 Nonlocking Screw Insertion in Distal Portion of Plate

Secure the plate to a distal fragment by inserting a 2.3 mm Nonlocking Cortical (Hex) Screw (CO-N23XX) through the medial-most center hole. Place the 2.0 mm Locking Drill Guide 4 mm–32 mm (80-0249) into the center hole and turn clockwise so that the guide fully threads into the plate. Insert the 2.0 mm Quick Release Drill (80-0318) and advance to the desired depth. Determine drill depth by referencing where the laser mark on the drill aligns with the measurement on the 2.0 mm locking drill guide. Remove the drill guide and use the 1.5 mm Hex Driver Tip (HPC-0015) with the Cruciform Driver Handle (MS-2210) to advance the 2.3 mm nonlocking screw until the screw head fully engages the plate.



Figure 9



Tray Components

Superior Distal Clavicle Plates

1	16-Hole Distal Clavicle Plate, 2.3 mm, Right	70-0123
2	13-Hole Distal Clavicle Plate, 2.3 mm, Right	70-0125
3	12-Hole Distal Clavicle Plate, 3.5 mm, Right	70-0111
4	9-Hole Distal Clavicle Plate, 3.5 mm, Right	70-0116
5	9-Hole Distal Clavicle Plate, 3.5 mm, Left	70-0117
6	12-Hole Distal Clavicle Plate, 3.5 mm, Left	70-0112
7	13-Hole Distal Clavicle Plate, 2.3 mm, Left	70-0126
8	16-Hole Distal Clavicle Plate, 2.3 mm, Left	70-0124
9	Low-Profile Clavicle J-Plate, 8-Hole, Left	70-0319
10	Low-Profile Clavicle J-Plate, 8-Hole, Right	70-0320

Optional Sterile Implants

16-Hole Distal Clavicle Plate, 3.5 mm, Left	7002-0416L-S
16-Hole Distal Clavicle Plate, 3.5 mm, Right	7002-0416R-S

Anterior Medial and Lateral Clavicle Plates

11	6-Hole Lateral Anterior Clavicle Plate	70-0122
12	8-Hole Lateral Anterior Clavicle Plate	70-0118
13	6-Hole Medial Anterior Clavicle Plate	70-0120
14	8-Hole Medial Anterior Clavicle Plate	70-0119
15	10-Hole Anterior Clavicle Plate	70-0121

Instrumentation

16	2.3 mm Screw Sleeve	MS-SS23
17	Cruciform Driver Handle	MS-2210
18	2.0 mm Quick Release Drill*	80-0318
19	2.0 mm Locking Drill Guide, 4 mm–32 mm	80-0249
20	Drill Guide for Distal Screws	MS-LDG23
21	1.5 mm Hex Driver Tip	HPC-0015
22	Targeting Guide, Distal Clavicle Plate, Left	80-0451
23	Targeting Guide, Distal Clavicle Plate, Right	80-0450


*Implants and screws are also available sterile-packed. Add an “-S” at end of product number for sterile product. For more details on sterile products, including pricing, contact our Business Services Department toll free at 888.627.9957.

Ordering Information [continued]


Sterile Plates			
Distal Clavicle Plate 3.5 mm 12-hole, Right	70-0111-S	Low-profile Clavicle Plate, 8-hole Large, Right	70-0289-S
Distal Clavicle Plate 3.5 mm 12-hole, Left	70-0112-S	Low-profile Clavicle Plate, 8-hole Medium, Left	70-0290-S
Distal Clavicle Plate 3.5 mm 9-hole, Right	70-0116-S	Low-profile Clavicle Plate, 8-hole Medium, Right	70-0291-S
Distal Clavicle Plate 3.5 mm 9-hole, Left	70-0117-S	Low-profile Clavicle Plate, 8-hole Small, Left	70-0292-S
8-hole Lateral Anterior Clavicle Plate	70-0118-S	Low-profile Clavicle Plate, 8-hole Small, Right	70-0293-S
8-hole Medial Anterior Clavicle Plate	70-0119-S	Low-profile Clavicle Plate, 10-hole Left	70-0294-S
6-hole Medial Anterior Clavicle Plate	70-0120-S	Low-profile Clavicle Plate, 10-hole Right	70-0295-S
10-hole Anterior Clavicle Plate	70-0121-S	Narrow-profile Clavicle Plate, 6-hole Left	70-0296-S
6-hole Lateral Anterior Clavicle Plate	70-0122-S	Narrow-profile Clavicle Plate, 6-hole Right	70-0297-S
Distal Clavicle Plate 2.3 mm 16-hole, Right	70-0123-S	Narrow-profile Clavicle Plate, 8-hole Straight, Left	70-0298-S
Distal Clavicle Plate 2.3 mm 16-hole, Left	70-0124-S	Narrow-profile Clavicle Plate, 8-hole Straight, Right	70-0299-S
Distal Clavicle Plate 2.3 mm 13-hole, Right	70-0125-S	Narrow-profile Clavicle Plate, 8-hole Large, Left	70-0300-S
Distal Clavicle Plate 2.3 mm 13-hole, Left	70-0126-S	Narrow-profile Clavicle Plate, 8-hole Large, Right	70-0301-S
Low-profile Clavicle Plate, 8-hole Straight, Left	70-0286-S	Low-profile Clavicle J-Plate, 8-hole Left	70-0319-S
Low-profile Clavicle Plate, 8-hole Straight, Right	70-0287-S	Low-profile Clavicle J-Plate, 8-hole Right	70-0320-S
Low-profile Clavicle Plate, 8-hole Large, Left	70-0288-S		

Clavicle Plating Diagrams

Anterior Medial and Lateral Clavicle Plates



Low and Narrow-profile Superior Midshaft Clavicle Plates



Superior Distal Clavicle Plates

