

Item Catalogue

NUPCO CODE: 4232150045600 | SN: 1233

NUPCO ITEM DESCRIPTION

SCREW, VARIABLE ANGLED 2.7 MM VARIABLE ANGLED CLAVICLE PLATING ASSORTED SIZES

VENDOR PRODUCT DESCRIPTION


Clavicle Plating System 2.7 mm locking screws

SKU CODES INCLUDED FOR THIS NUPCO CODE

- COL-2080
- COL-2100
- COL-2120
- COL-2140
- COL-2160
- COL-2180
- COL-2200
- COL-2220
- COL-2240
- COL-2260
- COL-2280
- COL-2300
- COL-2320
- COL-2340
- COL-2360
- COL-2380
- COL-2400
- COL-2420
- COL-2440
- COL-2460
- COL-2480
- COL-2500
- COL-2550
- COL-2600
- COL-2650

Superior Distal Clavicle Plate Surgical Technique [continued]

Figure 8



6 Nonlocking Screw Insertion in Distal Portion of Plate

Secure the plate to a distal fragment by inserting a 2.3 mm Nonlocking Cortical (Hex) Screw (CO-N23XX) through the medial-most center hole. Place the 2.0 mm Locking Drill Guide 4 mm–32 mm (80-0249) into the center hole and turn clockwise so that the guide fully threads into the plate. Insert the 2.0 mm Quick Release Drill (80-0318) and advance to the desired depth. Determine drill depth by referencing where the laser mark on the drill aligns with the measurement on the 2.0 mm locking drill guide. Remove the drill guide and use the 1.5 mm Hex Driver Tip (HPC-0015) with the Cruciform Driver Handle (MS-2210) to advance the 2.3 mm nonlocking screw until the screw head fully engages the plate.



Figure 9

Ordering Information [continued]

Screws			
2.7 mm (Nonlocking) Cortical (Hex) Screws*		2.7 mm Locking Cortical (Hex) Screws*	
2.7 mm x 8 mm Cortical Screw	CO-2708	2.7 mm x 8 mm Locking Cortical Screw	COL-2080
2.7 mm x 10 mm Cortical Screw	CO-2710	2.7 mm x 10 mm Locking Cortical Screw	COL-2100
2.7 mm x 12 mm Cortical Screw	CO-2712	2.7 mm x 12 mm Locking Cortical Screw	COL-2120
2.7 mm x 14 mm Cortical Screw	CO-2714	2.7 mm x 14 mm Locking Cortical Screw	COL-2140
2.7 mm x 16 mm Cortical Screw	CO-2716	2.7 mm x 16 mm Locking Cortical Screw	COL-2160
2.7 mm x 18 mm Cortical Screw	CO-2718	2.7 mm x 18 mm Locking Cortical Screw	COL-2180
2.7 mm x 20 mm Cortical Screw	CO-2720	2.7 mm x 20 mm Locking Cortical Screw	COL-2200
2.7 mm x 22 mm Cortical Screw	CO-2722	2.7 mm x 22 mm Locking Cortical Screw	COL-2220
2.7 mm x 24 mm Cortical Screw	CO-2724	2.7 mm x 24 mm Locking Cortical Screw	COL-2240
2.7 mm x 26 mm Cortical Screw	CO-2726	2.7 mm x 26 mm Locking Cortical Screw	COL-2260
2.7 mm x 28 mm Cortical Screw	CO-2728	2.7 mm x 28 mm Locking Cortical Screw	COL-2280
2.7 mm x 30 mm Cortical Screw	CO-2730	2.7 mm x 30 mm Locking Cortical Screw	COL-2300
2.7 mm x 32 mm Cortical Screw	CO-2732	2.7 mm x 32 mm Locking Cortical Screw	COL-2320
2.7 mm x 34 mm Cortical Screw	CO-2734	2.7 mm x 34 mm Locking Cortical Screw	COL-2340
2.7 mm x 36 mm Cortical Screw	CO-2736	2.7 mm x 36 mm Locking Cortical Screw	COL-2360
2.7 mm x 38 mm Cortical Screw	CO-2738	2.7 mm x 38 mm Locking Cortical Screw	COL-2380
2.7 mm x 40 mm Cortical Screw	CO-2740	2.7 mm x 40 mm Locking Cortical Screw	COL-2400
2.7 mm x 45 mm Cortical Screw	CO-2745	2.7 mm x 45 mm Locking Cortical Screw	COL-2450
2.7 mm x 50 mm Cortical Screw	CO-2750	2.7 mm x 50 mm Locking Cortical Screw	COL-2500
2.7 mm x 55 mm Cortical Screw	CO-2755	2.7 mm x 55 mm Locking Cortical Screw	COL-2550
2.7 mm x 60 mm Cortical Screw	CO-2760	2.7 mm x 60 mm Locking Cortical Screw	COL-2600
2.7 mm x 65 mm Cortical Screw	CO-2765	2.7 mm x 65 mm Locking Cortical Screw	COL-2650

*Implants and screws are also available sterile-packed. Add an “-S” at end of product number for sterile product. For more details on sterile products, including pricing, contact our Business Services Department toll free at 888.627.9957.

