

Item Catalogue

NUPCO CODE: 4232191208200 | SN: 1237

NUPCO ITEM DESCRIPTION

PLATE, VARIABLE ANGLED LCP CLAVICLE PLATES 2.7 MEDIAL ASSORTED SIZES, SHAPES, AND TYPES, AS REQUESTED BY END USER

VENDOR PRODUCT DESCRIPTION


Clavicle Plating System superior midshaft clavicle plates

SKU CODES INCLUDED FOR THIS NUPCO CODE

- 70-0286
- 70-0287
- 70-0288
- 70-0289
- 70-0290
- 70-0291
- 70-0292
- 70-0293
- 70-0294
- 70-0295
- 70-0296
- 70-0297
- 70-0298
- 70-0299
- 70-0300
- 70-0301

Superior Distal Clavicle Plate Surgical Technique [continued]

Figure 8




6 Nonlocking Screw Insertion in Distal Portion of Plate

Secure the plate to a distal fragment by inserting a 2.3 mm Nonlocking Cortical (Hex) Screw (CO-N23XX) through the medial-most center hole. Place the 2.0 mm Locking Drill Guide 4 mm–32 mm (80-0249) into the center hole and turn clockwise so that the guide fully threads into the plate. Insert the 2.0 mm Quick Release Drill (80-0318) and advance to the desired depth. Determine drill depth by referencing where the laser mark on the drill aligns with the measurement on the 2.0 mm locking drill guide. Remove the drill guide and use the 1.5 mm Hex Driver Tip (HPC-0015) with the Cruciform Driver Handle (MS-2210) to advance the 2.3 mm nonlocking screw until the screw head fully engages the plate.



Figure 9

Ordering Information



Tray Components

Superior Midshaft Clavicle Plates


1	Low-profile 10-hole, Right	70-0295	13	Narrow-profile, 6-hole, Left	70-0296
2	Low-profile 8-hole, Straight, Right	70-0287	14	Narrow-profile, 8-hole, Straight, Left	70-0298
3	Low-profile 8-hole, Large, Right	70-0289	15	Narrow-profile, 8-hole, Large, Left	70-0300
4	Low-profile 8-hole, Medium, Right	70-0291	16	Low-profile 10-hole, Left	70-0294
5	Low-profile 8-hole, Small, Right	70-0293	Instrumentation		
6	Narrow-profile, 6-hole, Right	70-0297	17	Clavicle Retractor	PL-CL03
7	Narrow-profile, 8-hole, Straight, Right	70-0299	18	CO/CA Countersink	PL-2080
8	Narrow-profile, 8-hole, Large, Right	70-0301	19	2.5 mm Solid, Quick Release, Driver Tip	HT-2502
9	Low-profile 8-hole, Straight, Left	70-0286	20	2.5 mm Flexible Hex Driver	80-0302
10	Low-profile 8-hole, Large, Left	70-0288	21	2.8 mm/3.5 mm Lag Guide	MS-DS2835
11	Low-profile 8-hole, Medium, Left	70-0290	22	Sharp Hook	PL-CL06
12	Low-profile 8-hole, Small, Left	70-0292	23	3.5 mm Tap Sleeve Assembly	PL-2190

Ordering Information [continued]


Sterile Plates			
Distal Clavicle Plate 3.5 mm 12-hole, Right	70-0111-S	Low-profile Clavicle Plate, 8-hole Large, Right	70-0289-S
Distal Clavicle Plate 3.5 mm 12-hole, Left	70-0112-S	Low-profile Clavicle Plate, 8-hole Medium, Left	70-0290-S
Distal Clavicle Plate 3.5 mm 9-hole, Right	70-0116-S	Low-profile Clavicle Plate, 8-hole Medium, Right	70-0291-S
Distal Clavicle Plate 3.5 mm 9-hole, Left	70-0117-S	Low-profile Clavicle Plate, 8-hole Small, Left	70-0292-S
8-hole Lateral Anterior Clavicle Plate	70-0118-S	Low-profile Clavicle Plate, 8-hole Small, Right	70-0293-S
8-hole Medial Anterior Clavicle Plate	70-0119-S	Low-profile Clavicle Plate, 10-hole Left	70-0294-S
6-hole Medial Anterior Clavicle Plate	70-0120-S	Low-profile Clavicle Plate, 10-hole Right	70-0295-S
10-hole Anterior Clavicle Plate	70-0121-S	Narrow-profile Clavicle Plate, 6-hole Left	70-0296-S
6-hole Lateral Anterior Clavicle Plate	70-0122-S	Narrow-profile Clavicle Plate, 6-hole Right	70-0297-S
Distal Clavicle Plate 2.3 mm 16-hole, Right	70-0123-S	Narrow-profile Clavicle Plate, 8-hole Straight, Left	70-0298-S
Distal Clavicle Plate 2.3 mm 16-hole, Left	70-0124-S	Narrow-profile Clavicle Plate, 8-hole Straight, Right	70-0299-S
Distal Clavicle Plate 2.3 mm 13-hole, Right	70-0125-S	Narrow-profile Clavicle Plate, 8-hole Large, Left	70-0300-S
Distal Clavicle Plate 2.3 mm 13-hole, Left	70-0126-S	Narrow-profile Clavicle Plate, 8-hole Large, Right	70-0301-S
Low-profile Clavicle Plate, 8-hole Straight, Left	70-0286-S	Low-profile Clavicle J-Plate, 8-hole Left	70-0319-S
Low-profile Clavicle Plate, 8-hole Straight, Right	70-0287-S	Low-profile Clavicle J-Plate, 8-hole Right	70-0320-S
Low-profile Clavicle Plate, 8-hole Large, Left	70-0288-S		

Clavicle Plating Diagrams

Anterior Medial and Lateral Clavicle Plates



Low and Narrow-profile Superior Midshaft Clavicle Plates



Superior Distal Clavicle Plates

