

Item Catalogue

NUPCO CODE: 4232191208300 | SN: 1238

NUPCO ITEM DESCRIPTION

PLATE, VARIABLE ANGLED LCP CLAVICLE HOOK PLATES 2.7 ASSORTED SIZES, SHAPES, AND TYPES, AS REQUESTED BY END USER

VENDOR PRODUCT DESCRIPTION

Clavicle Hook Plate implants

SKU CODES INCLUDED FOR THIS NUPCO CODE

7014-0712L
7014-0712R
7014-0716L
7014-0716R
7014-0720L
7014-0720R
7014-0612L
7014-0612R
7014-0616L
7014-0616R
7014-0620L
7014-0620R
7014-0512L
7014-0512R
7014-0516L
7014-0516R
7014-0520L
7014-0520R

Ordering Information

Tray Components

Clavicle Hook Plate Implants

1	Clavicle Hook Plate, 7-Hole 12 mm, Left	7014-0712L	10	Clavicle Hook Plate, 7-Hole 12 mm, Right	7014-0712R
2	Clavicle Hook Plate, 6-Hole 12 mm, Left	7014-0612L	11	Clavicle Hook Plate, 6-Hole 12 mm, Right	7014-0612R
3	Clavicle Hook Plate, 5-Hole 12 mm, Left	7014-0512L*	12	Clavicle Hook Plate, 5-Hole 12 mm, Right	7014-0512R*
4	Clavicle Hook Plate, 7-Hole 16 mm, Left	7014-0716L	13	Clavicle Hook Plate, 7-Hole 16 mm, Right	7014-0716R
5	Clavicle Hook Plate, 6-Hole 16 mm, Left	7014-0616L	14	Clavicle Hook Plate, 6-Hole 16 mm, Right	7014-0616R
6	Clavicle Hook Plate, 5-Hole 16 mm, Left	7014-0516L*	15	Clavicle Hook Plate, 5-Hole 16 mm, Right	7014-0516R*
7	Clavicle Hook Plate, 7-Hole 20 mm, Left	7014-0720L	16	Clavicle Hook Plate, 7-Hole 20 mm, Right	7014-0720R
8	Clavicle Hook Plate, 6-Hole 20 mm, Left	7014-0620L	17	Clavicle Hook Plate, 6-Hole 20 mm, Right	7014-0620R
9	Clavicle Hook Plate, 5-Hole 20 mm, Left	7014-0520L*	18	Clavicle Hook Plate, 5-Hole 20 mm, Right	7014-0520R*

Instrumentation

19	Clavicle Hook Plate Length Sizer	80-4007
20	Clavicle Hook Plate Depth Sizer, Left	80-4011
21	Clavicle Hook Plate Depth Sizer, Right	80-4012
22	Clavicle Hook Plate Depth Sizer, Hook	80-4010
23	3.0 mm Easyout, QR	80-0601
24	2.5 mm Easyout, QR	80-0600

Trays

25	Clavicle Hook Plate Tray	80-4020
26	Clavicle Hook Plate Essential Tray**	80-4019
27	Clavicle Hook Plate Tray Lid**	80-4021

*5-hole Clavicle Hook Plates are optional

** Not shown

Optional

Sterile Implants

Clavicle Hook Plate, 9-Hole 12 mm, Left	7014-0912L-S	Clavicle Hook Plate, 9-Hole 20 mm, Left	7014-0920L-S
Clavicle Hook Plate, 9-Hole 12 mm, Right	7014-0912R-S	Clavicle Hook Plate, 9-Hole 20 mm, Right	7014-0920R-S
Clavicle Hook Plate, 9-Hole 16 mm, Left	7014-0916L-S		
Clavicle Hook Plate, 9-Hole 16 mm, Right	7014-0916R-S		

Ordering Information [continued]

Sterile Clavicle Hook Plate Implants

Clavicle Hook Plate, 7-Hole 12 mm, Left	7014-0712L-S	Clavicle Hook Plate, 7-Hole 12 mm, Right	7014-0712R-S
Clavicle Hook Plate, 6-Hole 12 mm, Left	7014-0612L-S	Clavicle Hook Plate, 6-Hole 12 mm, Right	7014-0612R-S
Clavicle Hook Plate, 5-Hole 12 mm, Left	7014-0512L-S	Clavicle Hook Plate, 5-Hole 12 mm, Right	7014-0512R-S
Clavicle Hook Plate, 7-Hole 16 mm, Left	7014-0716L-S	Clavicle Hook Plate, 7-Hole 16 mm, Right	7014-0716R-S
Clavicle Hook Plate, 6-Hole 16 mm, Left	7014-0616L-S	Clavicle Hook Plate, 6-Hole 16 mm, Right	7014-0616R-S
Clavicle Hook Plate, 5-Hole 16 mm, Left	7014-0516L-S	Clavicle Hook Plate, 5-Hole 16 mm, Right	7014-0516R-S
Clavicle Hook Plate, 7-Hole 20 mm, Left	7014-0720L-S	Clavicle Hook Plate, 7-Hole 20 mm, Right	7014-0720R-S
Clavicle Hook Plate, 6-Hole 20 mm, Left	7014-0620L-S	Clavicle Hook Plate, 6-Hole 20 mm, Right	7014-0620R-S
Clavicle Hook Plate, 5-Hole 20 mm, Left	7014-0520L-S	Clavicle Hook Plate, 5-Hole 20 mm, Right	7014-0520R-S

Tray Components

Instrumentation

Plate Bender	PL-2040	2.8 mm Quick Release Drill	80-0387
Plate Bender, Large	PL-2045	.059" x 5" ST Guide Wire*	WS-1505ST
Depth Gauge 6–65 mm	80-0623	T15 Stick Fit Hexalobe Driver	80-0760
2.8 mm/3.5 mm Thin Drill Guide	PL-2196	Large Cannulated Quick Release Driver Handle	MS-3200
2.8 mm Hexalobe Locking Drill Guide 6–65 mm	80-0668	6 mm–70 mm Depth Gauge, 2 mm Increments	MS-9022
2.5 mm Quick Release Hex Driver	HPC-0025	Reduction Forceps w/Serrated Jaw	PL-CL04
3.5 mm Cortical Screw Bone Tap	MS-LTT35	3.5 mm Locking Drill Guide	MS-LDG35
Plate Tack	PL-PTACK	3.5 mm Screw Driver Sleeve	MS-SS35
3.5 mm x 5" Quick Release Drill	MS-DC35	Clavicle Retractor	PL-CL03
2.8 mm x 5" Quick Release Drill	MS-DC28		


*Also used as a K-wire

References

1. Kashii M, Inui H, Yamamoto K. Surgical treatment of distal clavicle fractures using the clavicular hook plate. *Clin Orthop Relat Res.* 2006;447:158-164.

Superior Distal Clavicle Plate Surgical Technique [continued]

Figure 8



6 Nonlocking Screw Insertion in Distal Portion of Plate

Secure the plate to a distal fragment by inserting a 2.3 mm Nonlocking Cortical (Hex) Screw (CO-N23XX) through the medial-most center hole. Place the 2.0 mm Locking Drill Guide 4 mm–32 mm (80-0249) into the center hole and turn clockwise so that the guide fully threads into the plate. Insert the 2.0 mm Quick Release Drill (80-0318) and advance to the desired depth. Determine drill depth by referencing where the laser mark on the drill aligns with the measurement on the 2.0 mm locking drill guide. Remove the drill guide and use the 1.5 mm Hex Driver Tip (HPC-0015) with the Cruciform Driver Handle (MS-2210) to advance the 2.3 mm nonlocking screw until the screw head fully engages the plate.



Figure 9