

Item Catalogue

NUPCO CODE: 4229520908500 | SN: 1484

NUPCO ITEM DESCRIPTION

CANNULATED DRILL BIT, ASSORTED SIZES AND LENGTHS FROM 2MM TO 3.0 MM

VENDOR PRODUCT DESCRIPTION

ExtremiFix Mini/Small cannulated drills 2.3/2.7 mm

SKU CODES INCLUDED FOR THIS NUPCO CODE

- 316-0529
- 316-0531
- 316-0532
- 316-0537
- 316-0538

Instrument Overview [continued]

2.0/2.4 mm Mini Instruments



**Ø1.7 mm Cannulated Drill,
Short, Quick Release**
(316-0513)



**Ø1.7 mm Cannulated Drill,
Long, Quick Release**
(316-0514)



**Ø2.0/2.4 mm Cannulated
Countersink, Quick Release**
(316-0511)



**Ø2.0/2.4 mm Proximal Cortex
Drill, Quick Release**
(316-0516)



**Ø1.7 mm Drill/Countersink
Combo, Short, Quick Release**
(316-0527)



**Ø1.7 mm Drill/Countersink
Combo, Long, Quick Release**
(316-0528)



**Ø2.0/2.4 mm T7 Driver,
Cannulated, Short,
Quick Release**
(316-0544)




**Ø2.0/2.4 mm T7 Driver,
Cannulated, Long,
Quick Release**
(316-0545)


Surgical Technique Overview

Cannulated Screw Surgical Technique

Reduction and
Provisional Fixation



Insert K-wire



Measure
Screw Length



Instrument Overview [continued]

3.0/4.0 mm Small Instruments



Ø3.0/4.0 mm Cannulated Countersink, Quick Release
(316-0501)



Ø3.0/4.0 mm T10 Driver, Cannulated, Short, Quick Release
(316-0554)



Ø3.0/4.0 mm T10 Driver, Cannulated, Long, Quick Release
(316-0555)

3.0 mm Small Instruments



Ø2.3 mm Cannulated Drill, Quick Release
(316-0532)



Ø3.2 mm Proximal Cortex Drill, Quick Release
(316-0533)



Ø2.3 mm Drill/Countersink Combo, Short, Quick Release
(316-0537)



Ø2.3 mm Drill/Countersink Combo, Long, Quick Release
(316-0538)

4.0 mm Small Instruments



Ø2.7 mm Cannulated Drill, Quick Release
(316-0534)



Ø3.7 mm Proximal Cortex Drill, Quick Release
(316-0543)




Ø2.7 mm Drill/Countersink Combo, Short, Quick Release
(316-0547)



Ø2.7 mm Drill/Countersink Combo, Long, Quick Release
(316-0548)

Predrill (Recommended)



Countersink or Cortex Drill (Recommended)



Insert Screw



Screw Removal (Optional)



Foot & Ankle Applications

Headed and Headless Cannulated Screws [continued]




Figure 3

3 Measure Screw Length

Slide the Cannulated Depth Gauge, Double Sided (316-0036) over the K-wire until the tip of the depth gauge contacts the surface of the bone (Figure 3). Make sure to use the side of the depth gauge that corresponds to the length of guide wire used.

Note: After screw length has been recorded and when appropriate, the guide wire may be advanced to maintain position during pre-drilling.




Figure 4

4 Predrill (Recommended)

The ExtremiFix Mini & Small cannulated screws are self-drilling and self-tapping.

Note: In cases of dense cortical bone, pre-drilling with the appropriately sized cannulated drill is recommended.

Drills for 2.0/2.4 mm Screws

316-0513	Ø1.7 mm Cannulated Drill, Short, Quick Release
316-0514	Ø1.7 mm Cannulated Drill, Long, Quick Release

Drill for 3.0/4.0 mm Screws

316-0532	Ø2.3 mm Cannulated Drill, Quick Release
-----------------	---

Use the selected Cannulated Drill (316-05XX) through the soft tissue protector to drill a pilot hole (Figure 4).

Caution: To maintain K-wire position, drilling should not go beyond the tip of the wire.



Cannulated Depth Gauge, Double Sided (316-0036)



K-wire (316-0XXX)



Ø1.7 mm Cannulated Drill, Short, Quick Release (316-0513)



Ø1.7 mm Cannulated Drill, Long, Quick Release (316-0514)



Ø2.3 mm Cannulated Drill, Quick Release (316-0532)

Ordering Information [continued]

Tray Components — 3.0 mm Driver and Drill Module

1	Ø3.0 mm Screw & Instrument Module, ExtremiFix Small	316-1530	2	Ø2.3 mm Cannulated Drill, Quick Release	316-0532
3	Ø3.0/4.0 mm Cannulated Countersink, Quick Release	316-0501	4	Ø3.2 mm Proximal Cortex Drill, Quick Release	316-0533
5	Ø2.3 mm Drill/Countersink Combo, Short, Quick Release	316-0537	6	Ø2.3 mm Drill/Countersink Combo, Long, Quick Release	316-0538
7	Ø3.0/4.0 mm T10 Driver, Cannulated, Short, Quick Release	316-0554	8	Ø3.0/4.0 mm T10 Driver, Cannulated, Long, Quick Release	316-0555


Tray Components — 3.0 mm Screw Module

9	Ø3.0 mm Screw & Instrument Module, ExtremiFix Small	316-1530	10	Cannulated Depth Gauge, Double sided 4"/6"	316-0036
11	Ø.045 x 6" K-wire	316-0004	12	Ø.045 x 6" Partially Threaded K-wire	316-0006
13	Ø.045 x 6" Double Trocar K-wire	316-0122	14	Ø.045 x 4" K-wire	316-0003
15	Ø.045 x 4" Double Trocar K-wire	316-0123	16	Ø3.0 mm Washer	316-0225

Tray Components — 3.0 mm Screw Module

17 3.0 mm Headed Cannulated Screws		18 3.0 mm Headless Cannulated Screws	
3.0 x 10 mm Cannulated Headed Screw	318-3010	3.0 x 10 mm Cannulated Headless Screw	315-3010
3.0 x 12 mm Cannulated Headed Screw	318-3012	3.0 x 12 mm Cannulated Headless Screw	315-3012
3.0 x 14 mm Cannulated Headed Screw	318-3014	3.0 x 14 mm Cannulated Headless Screw	315-3014
3.0 x 16 mm Cannulated Headed Screw	318-3016	3.0 x 16 mm Cannulated Headless Screw	315-3016
3.0 x 18 mm Cannulated Headed Screw	318-3018	3.0 x 18 mm Cannulated Headless Screw	315-3018
3.0 x 20 mm Cannulated Headed Screw	318-3020	3.0 x 20 mm Cannulated Headless Screw	315-3020
3.0 x 22 mm Cannulated Headed Screw	318-3022	3.0 x 22 mm Cannulated Headless Screw	315-3022
3.0 x 24 mm Cannulated Headed Screw	318-3024	3.0 x 24 mm Cannulated Headless Screw	315-3024
3.0 x 26 mm Cannulated Headed Screw	318-3026	3.0 x 26 mm Cannulated Headless Screw	315-3026
3.0 x 28 mm Cannulated Headed Screw	318-3028	3.0 x 28 mm Cannulated Headless Screw	315-3028
3.0 x 30 mm Cannulated Headed Screw	318-3030	3.0 x 30 mm Cannulated Headless Screw	315-3030
3.0 x 32 mm Cannulated Headed Screw	318-3032	3.0 x 32 mm Cannulated Headless Screw	315-3032
3.0 x 34 mm Cannulated Headed Screw	318-3034	3.0 x 34 mm Cannulated Headless Screw	315-3034
3.0 x 36 mm Cannulated Headed Screw	318-3036	3.0 x 36 mm Cannulated Headless Screw	315-3036
3.0 x 38 mm Cannulated Headed Screw	318-3038	3.0 x 38 mm Cannulated Headless Screw	315-3038
3.0 x 40 mm Cannulated Headed Screw	318-3040	3.0 x 40 mm Cannulated Headless Screw	315-3040

Ordering Information [continued]



Instrument Overview

Flip-up Caddy Instruments



1.3 mm Threaded Guide Wire, 150 mm
(80-2038)



1.3 mm Smooth Guide Wire, 150 mm
(80-2039)



1.3 mm Wire Depth Gauge
(80-3960)



2.7 mm Cannulated Drill, Quick Connect
(80-2075)



Cannula Handle
(80-3961)



1.3 mm Wire Cannula
(80-3962)



2.7 mm Drill Cannula
(80-3963)



4.0 mm Screw Cannula
(80-3964)



2.5 mm Short Cannulated Hex Driver
(80-3956)



5.3 mm Cannulated Countersink
(80-2042)



4.0 mm Cannulated Tap
(80-2081)



1.25 mm Cleaning Stylet
(80-2246)



Large Screw Holding Forceps
(MS-45210)



2.5 mm Screw Driver Sleeve
(80-3957)



1.3/2.7 mm Double Drill Sleeve
(80-2079)



1.3 mm Parallel Guide Base
(80-3778)



1.3 mm Parallel Guide Drop-in Cannula
(80-3781)



Small Ratchet Handle with Quick Release Connection
(80-0398)



2.5 mm Solid Hex Driver, Quick Connect
(80-2074)



2.5 mm Easyout, Quick Release
(80-0600)

