

Item Catalogue

NUPCO CODE: 4232150020100 | SN: 685

NUPCO ITEM DESCRIPTION

SCREW, LOCKED HEAD 3.5MM OR EQUIVALENT FOR STAINLESS STEEL HOOK CLAVICLE PLATE

VENDOR PRODUCT DESCRIPTION

Acumed Clavicle Hook Plate 3.5 mm locking screw families

SKU CODES INCLUDED FOR THIS NUPCO CODE

- 30-0232
- 30-0233
- 30-0234
- 30-0235
- 30-0236
- 30-0237
- 30-0238
- 30-0239
- 30-0240
- 30-0241
- 30-0242
- 30-0243
- 30-0244
- 30-0245
- 30-0246
- 30-0247
- 30-0248
- 30-3642
- 30-3644
- 30-3646
- 30-3648
- 30-0250
- 30-0251
- 30-0252
- 30-0253

Ordering Information [continued]

Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264

3.5 mm (Nonlocking) Cortical (Hex) Screws

3.5 mm x 8 mm Cortical Screw	CO-3080
3.5 mm x 10 mm Cortical Screw	CO-3100
3.5 mm x 12 mm Cortical Screw	CO-3120
3.5 mm x 14 mm Cortical Screw	CO-3140
3.5 mm x 16 mm Cortical Screw	CO-3160
3.5 mm x 18 mm Cortical Screw	CO-3180
3.5 mm x 20 mm Cortical Screw	CO-3200
3.5 mm x 22 mm Cortical Screw	CO-3220
3.5 mm x 24 mm Cortical Screw	CO-3240
3.5 mm x 26 mm Cortical Screw	CO-3260

3.5 mm Locking Cortical (Hex) Screws

3.5 mm x 8 mm Locking Cortical Screw	COL-3080
3.5 mm x 10 mm Locking Cortical Screw	COL-3100
3.5 mm x 12 mm Locking Cortical Screw	COL-3120
3.5 mm x 14 mm Locking Cortical Screw	COL-3140
3.5 mm x 16 mm Locking Cortical Screw	COL-3160
3.5 mm x 18 mm Locking Cortical Screw	COL-3180
3.5 mm x 20 mm Locking Cortical Screw	COL-3200
3.5 mm x 22 mm Locking Cortical Screw	COL-3220
3.5 mm x 24 mm Locking Cortical Screw	COL-3240
3.5 mm x 26 mm Locking Cortical Screw	COL-3260

Ordering Information [continued]

Sterile Screws

Sterile 3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232-S
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233-S
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234-S
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235-S
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236-S
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237-S
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238-S
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239-S
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240-S
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241-S

Sterile 3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255-S
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256-S
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257-S
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258-S
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259-S
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260-S
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261-S
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262-S
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263-S
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264-S

Sterile 3.5 mm (Nonlocking) Cortical (Hex) Screws

3.5 mm x 8 mm Cortical Screw	CO-3080-S
3.5 mm x 10 mm Cortical Screw	CO-3100-S
3.5 mm x 12 mm Cortical Screw	CO-3120-S
3.5 mm x 14 mm Cortical Screw	CO-3140-S
3.5 mm x 16 mm Cortical Screw	CO-3160-S
3.5 mm x 18 mm Cortical Screw	CO-3180-S
3.5 mm x 20 mm Cortical Screw	CO-3200-S
3.5 mm x 22 mm Cortical Screw	CO-3220-S
3.5 mm x 24 mm Cortical Screw	CO-3240-S
3.5 mm x 26 mm Cortical Screw	CO-3260-S

Sterile 3.5 mm Locking Cortical (Hex) Screws

3.5 mm x 8 mm Locking Cortical Screw	COL-3080-S
3.5 mm x 10 mm Locking Cortical Screw	COL-3100-S
3.5 mm x 12 mm Locking Cortical Screw	COL-3120-S
3.5 mm x 14 mm Locking Cortical Screw	COL-3140-S
3.5 mm x 16 mm Locking Cortical Screw	COL-3160-S
3.5 mm x 18 mm Locking Cortical Screw	COL-3180-S
3.5 mm x 20 mm Locking Cortical Screw	COL-3200-S
3.5 mm x 22 mm Locking Cortical Screw	COL-3220-S
3.5 mm x 24 mm Locking Cortical Screw	COL-3240-S
3.5 mm x 26 mm Locking Cortical Screw	COL-3260-S

Note: To learn more about the full line of Acumed innovative surgical solutions, please contact your authorized Acumed distributor, call 888.627.9957, or visit www.acumed.net.

Superior Distal Clavicle Plate Surgical Technique [continued]

4 Plate Placement


Once the plate's ideal positioning has been selected, it is provisionally stabilized to the clavicle with .045" x 6" ST Guide Wire (WS-1106ST) or .059" x 5" ST Guide Wire (WS-1505ST) K-wires. Under radiographic evaluation, the most lateral K-wire hole of each Superior Distal Clavicle Plate (70-0XXX) affords the opportunity to verify that the placement of the screws will not protrude into the AC joint by inserting a K-wire to confirm plate placement.

Note: The Reduction Forceps With Serrated Jaw (PL-CL04) should be used for plate placement and are not designed to be used to reduce the plate to the bone or to hold the plate while attempting to bend or contour it to match the patient's anatomy. Plate Benders (PL-204X) are available in the event that plate contouring is required to achieve the desired fit to the clavicle.

Caution: If bending of the plate is necessary, please observe the following:

- ▶ Do not bend plates more than 30 degrees
- ▶ Bend radii should be greater than 1 inch
- ▶ Do not bend more than once
- ▶ Avoid bending across locking holes

Figure 5



Ordering Information [continued]

Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264

Hexalobe Short Screw Caddies

3.5 mm Hexalobe Short Screw Caddy	80-0843
3.5 mm Hexalobe Short Screw Caddy Lid	80-0856
3.0 mm Hexalobe Short Screw Caddy	80-1066
3.0 mm Hexalobe Short Screw Caddy Lid	80-1067

3.0 mm Locking Hexalobe Screws

3.0 mm x 8 mm Locking Hexalobe Screw	30-0278
3.0 mm x 10 mm Locking Hexalobe Screw	30-0279
3.0 mm x 12 mm Locking Hexalobe Screw	30-0280
3.0 mm x 14 mm Locking Hexalobe Screw	30-0281
3.0 mm x 16 mm Locking Hexalobe Screw	30-0282
3.0 mm x 18 mm Locking Hexalobe Screw	30-0283
3.0 mm x 20 mm Locking Hexalobe Screw	30-0284
3.0 mm x 22 mm Locking Hexalobe Screw	30-0285
3.0 mm x 24 mm Locking Hexalobe Screw	30-0286
3.0 mm x 26 mm Locking Hexalobe Screw	30-0287

3.0 mm Nonlocking Hexalobe Screws

3.0 mm x 8 mm Nonlocking Hexalobe Screw	30-0301
3.0 mm x 10 mm Nonlocking Hexalobe Screw	30-0302
3.0 mm x 12 mm Nonlocking Hexalobe Screw	30-0303
3.0 mm x 14 mm Nonlocking Hexalobe Screw	30-0304
3.0 mm x 16 mm Nonlocking Hexalobe Screw	30-0305
3.0 mm x 18 mm Nonlocking Hexalobe Screw	30-0306
3.0 mm x 20 mm Nonlocking Hexalobe Screw	30-0307
3.0 mm x 22 mm Nonlocking Hexalobe Screw	30-0308
3.0 mm x 24 mm Nonlocking Hexalobe Screw	30-0309
3.0 mm x 26 mm Nonlocking Hexalobe Screw	30-0310

Ordering Information [continued]

Sterile Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232-S
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233-S
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234-S
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235-S
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236-S
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237-S
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238-S
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239-S
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240-S
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241-S

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255-S
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256-S
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257-S
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258-S
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259-S
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260-S
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261-S
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262-S
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263-S
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264-S

3.0 mm Locking Hexalobe Screws

3.0 mm x 8 mm Locking Hexalobe Screw	30-0278-S
3.0 mm x 10 mm Locking Hexalobe Screw	30-0279-S
3.0 mm x 12 mm Locking Hexalobe Screw	30-0280-S
3.0 mm x 14 mm Locking Hexalobe Screw	30-0281-S
3.0 mm x 16 mm Locking Hexalobe Screw	30-0282-S
3.0 mm x 18 mm Locking Hexalobe Screw	30-0283-S
3.0 mm x 20 mm Locking Hexalobe Screw	30-0284-S
3.0 mm x 22 mm Locking Hexalobe Screw	30-0285-S
3.0 mm x 24 mm Locking Hexalobe Screw	30-0286-S
3.0 mm x 26 mm Locking Hexalobe Screw	30-0287-S

3.0 mm Nonlocking Hexalobe Screws

3.0 mm x 8 mm Nonlocking Hexalobe Screw	30-0301-S
3.0 mm x 10 mm Nonlocking Hexalobe Screw	30-0302-S
3.0 mm x 12 mm Nonlocking Hexalobe Screw	30-0303-S
3.0 mm x 14 mm Nonlocking Hexalobe Screw	30-0304-S
3.0 mm x 16 mm Nonlocking Hexalobe Screw	30-0305-S
3.0 mm x 18 mm Nonlocking Hexalobe Screw	30-0306-S
3.0 mm x 20 mm Nonlocking Hexalobe Screw	30-0307-S
3.0 mm x 22 mm Nonlocking Hexalobe Screw	30-0308-S
3.0 mm x 24 mm Nonlocking Hexalobe Screw	30-0309-S
3.0 mm x 26 mm Nonlocking Hexalobe Screw	30-0310-S

Ordering Information [continued]

Screws

4 3.5 mm Nonlocking Hexalobe Screws*

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264
3.5 mm x 28 mm Nonlocking Hexalobe Screw	30-0265
3.5 mm x 30 mm Nonlocking Hexalobe Screw	30-0266
3.5 mm x 32 mm Nonlocking Hexalobe Screw	30-0267
3.5 mm x 34 mm Nonlocking Hexalobe Screw	30-0268
3.5 mm x 36 mm Nonlocking Hexalobe Screw	30-0269
3.5 mm x 38 mm Nonlocking Hexalobe Screw	30-0270
3.5 mm x 40 mm Nonlocking Hexalobe Screw	30-0271
3.5 mm x 42 mm Nonlocking Hexalobe Screw	30-3542
3.5 mm x 44 mm Nonlocking Hexalobe Screw	30-3544
3.5 mm x 46 mm Nonlocking Hexalobe Screw	30-3546
3.5 mm x 48 mm Nonlocking Hexalobe Screw	30-3548
3.5 mm x 50 mm Nonlocking Hexalobe Screw	30-0273
3.5 mm x 55 mm Nonlocking Hexalobe Screw	30-0274
3.5 mm x 60 mm Nonlocking Hexalobe Screw	30-0275
3.5 mm x 65 mm Nonlocking Hexalobe Screw	30-0276

5 3.5 mm Locking Hexalobe Screws*

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241
3.5 mm x 28 mm Locking Hexalobe Screw	30-0242
3.5 mm x 30 mm Locking Hexalobe Screw	30-0243
3.5 mm x 32 mm Locking Hexalobe Screw	30-0244
3.5 mm x 34 mm Locking Hexalobe Screw	30-0245
3.5 mm x 36 mm Locking Hexalobe Screw	30-0246
3.5 mm x 38 mm Locking Hexalobe Screw	30-0247
3.5 mm x 40 mm Locking Hexalobe Screw	30-0248
3.5 mm x 42 mm Locking Hexalobe Screw	30-3642
3.5 mm x 44 mm Locking Hexalobe Screw	30-3644
3.5 mm x 46 mm Locking Hexalobe Screw	30-3646
3.5 mm x 48 mm Locking Hexalobe Screw	30-3648
3.5 mm x 50 mm Locking Hexalobe Screw	30-0250
3.5 mm x 55 mm Locking Hexalobe Screw	30-0251
3.5 mm x 60 mm Locking Hexalobe Screw	30-0252
3.5 mm x 65 mm Locking Hexalobe Screw	30-0253