

Item Catalogue

NUPCO CODE: 4232150018500 | SN: 686

NUPCO ITEM DESCRIPTION

SCREW, CORTICAL 3.5MM OR EQUIVALENT FOR STAINLESS STEEL HOOK CLAVICLE PLATE, ASSORTED SIZES

VENDOR PRODUCT DESCRIPTION

Acumed Clavicle Hook Plate 3.5 mm cortical/nonlocking screw families

SKU CODES INCLUDED FOR THIS NUPCO CODE

- 30-0255
- 30-0256
- 30-0257
- 30-0258
- 30-0259
- 30-0260
- 30-0261
- 30-0262
- 30-0263
- 30-0264
- 30-0265
- 30-0266
- 30-0267
- 30-0268
- 30-0269
- 30-0270
- 30-0271
- 30-3542
- 30-3544
- 30-3546
- 30-3548
- 30-0273
- 30-0274
- 30-0275
- 30-0276

Ordering Information [continued]

Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264

3.5 mm (Nonlocking) Cortical (Hex) Screws

3.5 mm x 8 mm Cortical Screw	CO-3080
3.5 mm x 10 mm Cortical Screw	CO-3100
3.5 mm x 12 mm Cortical Screw	CO-3120
3.5 mm x 14 mm Cortical Screw	CO-3140
3.5 mm x 16 mm Cortical Screw	CO-3160
3.5 mm x 18 mm Cortical Screw	CO-3180
3.5 mm x 20 mm Cortical Screw	CO-3200
3.5 mm x 22 mm Cortical Screw	CO-3220
3.5 mm x 24 mm Cortical Screw	CO-3240
3.5 mm x 26 mm Cortical Screw	CO-3260

3.5 mm Locking Cortical (Hex) Screws

3.5 mm x 8 mm Locking Cortical Screw	COL-3080
3.5 mm x 10 mm Locking Cortical Screw	COL-3100
3.5 mm x 12 mm Locking Cortical Screw	COL-3120
3.5 mm x 14 mm Locking Cortical Screw	COL-3140
3.5 mm x 16 mm Locking Cortical Screw	COL-3160
3.5 mm x 18 mm Locking Cortical Screw	COL-3180
3.5 mm x 20 mm Locking Cortical Screw	COL-3200
3.5 mm x 22 mm Locking Cortical Screw	COL-3220
3.5 mm x 24 mm Locking Cortical Screw	COL-3240
3.5 mm x 26 mm Locking Cortical Screw	COL-3260

Ordering Information [continued]

Sterile Screws

Sterile 3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232-S
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233-S
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234-S
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235-S
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236-S
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237-S
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238-S
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239-S
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240-S
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241-S

Sterile 3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255-S
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256-S
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257-S
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258-S
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259-S
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260-S
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261-S
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262-S
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263-S
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264-S

Sterile 3.5 mm (Nonlocking) Cortical (Hex) Screws

3.5 mm x 8 mm Cortical Screw	CO-3080-S
3.5 mm x 10 mm Cortical Screw	CO-3100-S
3.5 mm x 12 mm Cortical Screw	CO-3120-S
3.5 mm x 14 mm Cortical Screw	CO-3140-S
3.5 mm x 16 mm Cortical Screw	CO-3160-S
3.5 mm x 18 mm Cortical Screw	CO-3180-S
3.5 mm x 20 mm Cortical Screw	CO-3200-S
3.5 mm x 22 mm Cortical Screw	CO-3220-S
3.5 mm x 24 mm Cortical Screw	CO-3240-S
3.5 mm x 26 mm Cortical Screw	CO-3260-S

Sterile 3.5 mm Locking Cortical (Hex) Screws

3.5 mm x 8 mm Locking Cortical Screw	COL-3080-S
3.5 mm x 10 mm Locking Cortical Screw	COL-3100-S
3.5 mm x 12 mm Locking Cortical Screw	COL-3120-S
3.5 mm x 14 mm Locking Cortical Screw	COL-3140-S
3.5 mm x 16 mm Locking Cortical Screw	COL-3160-S
3.5 mm x 18 mm Locking Cortical Screw	COL-3180-S
3.5 mm x 20 mm Locking Cortical Screw	COL-3200-S
3.5 mm x 22 mm Locking Cortical Screw	COL-3220-S
3.5 mm x 24 mm Locking Cortical Screw	COL-3240-S
3.5 mm x 26 mm Locking Cortical Screw	COL-3260-S

Note: To learn more about the full line of Acumed innovative surgical solutions, please contact your authorized Acumed distributor, call 888.627.9957, or visit www.acumed.net.

Superior Distal Clavicle Plate Surgical Technique [continued]

4 Plate Placement

Once the plate's ideal positioning has been selected, it is provisionally stabilized to the clavicle with .045" x 6" ST Guide Wire (WS-1106ST) or .059" x 5" ST Guide Wire (WS-1505ST) K-wires. Under radiographic evaluation, the most lateral K-wire hole of each Superior Distal Clavicle Plate (70-0XXX) affords the opportunity to verify that the placement of the screws will not protrude into the AC joint by inserting a K-wire to confirm plate placement.

Note: The Reduction Forceps With Serrated Jaw (PL-CL04) should be used for plate placement and are not designed to be used to reduce the plate to the bone or to hold the plate while attempting to bend or contour it to match the patient's anatomy. Plate Benders (PL-204X) are available in the event that plate contouring is required to achieve the desired fit to the clavicle.

Caution: If bending of the plate is necessary, please observe the following:

- ▶ Do not bend plates more than 30 degrees
- ▶ Bend radii should be greater than 1 inch
- ▶ Do not bend more than once
- ▶ Avoid bending across locking holes

Figure 5



Ordering Information [continued]

Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264

Hexalobe Short Screw Caddies

3.5 mm Hexalobe Short Screw Caddy	80-0843
3.5 mm Hexalobe Short Screw Caddy Lid	80-0856
3.0 mm Hexalobe Short Screw Caddy	80-1066
3.0 mm Hexalobe Short Screw Caddy Lid	80-1067

3.0 mm Locking Hexalobe Screws

3.0 mm x 8 mm Locking Hexalobe Screw	30-0278
3.0 mm x 10 mm Locking Hexalobe Screw	30-0279
3.0 mm x 12 mm Locking Hexalobe Screw	30-0280
3.0 mm x 14 mm Locking Hexalobe Screw	30-0281
3.0 mm x 16 mm Locking Hexalobe Screw	30-0282
3.0 mm x 18 mm Locking Hexalobe Screw	30-0283
3.0 mm x 20 mm Locking Hexalobe Screw	30-0284
3.0 mm x 22 mm Locking Hexalobe Screw	30-0285
3.0 mm x 24 mm Locking Hexalobe Screw	30-0286
3.0 mm x 26 mm Locking Hexalobe Screw	30-0287

3.0 mm Nonlocking Hexalobe Screws

3.0 mm x 8 mm Nonlocking Hexalobe Screw	30-0301
3.0 mm x 10 mm Nonlocking Hexalobe Screw	30-0302
3.0 mm x 12 mm Nonlocking Hexalobe Screw	30-0303
3.0 mm x 14 mm Nonlocking Hexalobe Screw	30-0304
3.0 mm x 16 mm Nonlocking Hexalobe Screw	30-0305
3.0 mm x 18 mm Nonlocking Hexalobe Screw	30-0306
3.0 mm x 20 mm Nonlocking Hexalobe Screw	30-0307
3.0 mm x 22 mm Nonlocking Hexalobe Screw	30-0308
3.0 mm x 24 mm Nonlocking Hexalobe Screw	30-0309
3.0 mm x 26 mm Nonlocking Hexalobe Screw	30-0310

Ordering Information [continued]

Sterile Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232-S
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233-S
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234-S
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235-S
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236-S
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237-S
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238-S
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239-S
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240-S
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241-S

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255-S
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256-S
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257-S
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258-S
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259-S
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260-S
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261-S
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262-S
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263-S
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264-S

3.0 mm Locking Hexalobe Screws

3.0 mm x 8 mm Locking Hexalobe Screw	30-0278-S
3.0 mm x 10 mm Locking Hexalobe Screw	30-0279-S
3.0 mm x 12 mm Locking Hexalobe Screw	30-0280-S
3.0 mm x 14 mm Locking Hexalobe Screw	30-0281-S
3.0 mm x 16 mm Locking Hexalobe Screw	30-0282-S
3.0 mm x 18 mm Locking Hexalobe Screw	30-0283-S
3.0 mm x 20 mm Locking Hexalobe Screw	30-0284-S
3.0 mm x 22 mm Locking Hexalobe Screw	30-0285-S
3.0 mm x 24 mm Locking Hexalobe Screw	30-0286-S
3.0 mm x 26 mm Locking Hexalobe Screw	30-0287-S

3.0 mm Nonlocking Hexalobe Screws

3.0 mm x 8 mm Nonlocking Hexalobe Screw	30-0301-S
3.0 mm x 10 mm Nonlocking Hexalobe Screw	30-0302-S
3.0 mm x 12 mm Nonlocking Hexalobe Screw	30-0303-S
3.0 mm x 14 mm Nonlocking Hexalobe Screw	30-0304-S
3.0 mm x 16 mm Nonlocking Hexalobe Screw	30-0305-S
3.0 mm x 18 mm Nonlocking Hexalobe Screw	30-0306-S
3.0 mm x 20 mm Nonlocking Hexalobe Screw	30-0307-S
3.0 mm x 22 mm Nonlocking Hexalobe Screw	30-0308-S
3.0 mm x 24 mm Nonlocking Hexalobe Screw	30-0309-S
3.0 mm x 26 mm Nonlocking Hexalobe Screw	30-0310-S

Variable Angle Screw Surgical Technique [continued]

4 Insert Variable Angle Screw

Note: Final tightening of the 2.7 mm and 3.5 mm variable angle hexalobe screws must be done manually and not under power. The Torque Limiting Quick Connect ensures a consistent insertion torque to provide a uniform screw-plate interface and may prevent overtightening of the screw. See below for the torque limit for each screw.

To insert a 2.7 mm variable angle hexalobe screw, assemble the 1.70 N•m Torque Limiting Quick Connect (80-2366) to the Handle for Torque Limiting Quick Connect (80-2368). Connect the T8 Stick Fit Hexalobe Driver (80-0759) to the Torque Limiting Quick Connect assembly.

To insert a 3.5 mm variable angle hexalobe screw, assemble the 2.26 N•m Torque Limiting Quick Connect (80-2367) to the Handle for Torque Limiting Quick Connect (80-2368). Connect the T15 Stick Fit Hexalobe Driver (80-0760) to the Torque Limiting Quick Connect assembly.

Advance the screw by hand until achieving an audible click and/or tactile feedback. Final tightening should be completed with the Torque Limiting Handle, which is designed to provide a secure lock between the plate and screw. Upon final seating, confirm proper screw placement and screw length under fluoroscopy.

Note: Do not use a Torque Limiting Quick Connect for screw removal.

Caution: Use the maximum number of screws based on the indication to reduce the risk of screw breakage during healing.




Figure 18

Screw	Torque Limit	Color Band
2.7 mm Variable Angle Hexalobe Screw (3013-27XXX)	1.70 N•m	Brown
3.5 mm Variable Angle Hexalobe Screw (3013-35XXX)	2.26 N•m	Black



Handle for Torque Limiting Quick Connect (80-2368)



T8 Stick Fit Hexalobe Driver (80-0759)



T15 Stick Fit Hexalobe Driver (80-0760)



2.7 mm Variable Angle Hexalobe Screw, 10–60 mm lengths (3013-27XXX)



3.5 mm Variable Angle Hexalobe Screw, 10–65 mm lengths (3013-35XXX)



1.70 N•m Torque Limiting Quick Connect (80-2366)



2.26 N•m Torque Limiting Quick Connect (80-2367)

Ordering Information [continued]

Screws

4 3.5 mm Nonlocking Hexalobe Screws*

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264
3.5 mm x 28 mm Nonlocking Hexalobe Screw	30-0265
3.5 mm x 30 mm Nonlocking Hexalobe Screw	30-0266
3.5 mm x 32 mm Nonlocking Hexalobe Screw	30-0267
3.5 mm x 34 mm Nonlocking Hexalobe Screw	30-0268
3.5 mm x 36 mm Nonlocking Hexalobe Screw	30-0269
3.5 mm x 38 mm Nonlocking Hexalobe Screw	30-0270
3.5 mm x 40 mm Nonlocking Hexalobe Screw	30-0271
3.5 mm x 42 mm Nonlocking Hexalobe Screw	30-3542
3.5 mm x 44 mm Nonlocking Hexalobe Screw	30-3544
3.5 mm x 46 mm Nonlocking Hexalobe Screw	30-3546
3.5 mm x 48 mm Nonlocking Hexalobe Screw	30-3548
3.5 mm x 50 mm Nonlocking Hexalobe Screw	30-0273
3.5 mm x 55 mm Nonlocking Hexalobe Screw	30-0274
3.5 mm x 60 mm Nonlocking Hexalobe Screw	30-0275
3.5 mm x 65 mm Nonlocking Hexalobe Screw	30-0276

5 3.5 mm Locking Hexalobe Screws*

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241
3.5 mm x 28 mm Locking Hexalobe Screw	30-0242
3.5 mm x 30 mm Locking Hexalobe Screw	30-0243
3.5 mm x 32 mm Locking Hexalobe Screw	30-0244
3.5 mm x 34 mm Locking Hexalobe Screw	30-0245
3.5 mm x 36 mm Locking Hexalobe Screw	30-0246
3.5 mm x 38 mm Locking Hexalobe Screw	30-0247
3.5 mm x 40 mm Locking Hexalobe Screw	30-0248
3.5 mm x 42 mm Locking Hexalobe Screw	30-3642
3.5 mm x 44 mm Locking Hexalobe Screw	30-3644
3.5 mm x 46 mm Locking Hexalobe Screw	30-3646
3.5 mm x 48 mm Locking Hexalobe Screw	30-3648
3.5 mm x 50 mm Locking Hexalobe Screw	30-0250
3.5 mm x 55 mm Locking Hexalobe Screw	30-0251
3.5 mm x 60 mm Locking Hexalobe Screw	30-0252
3.5 mm x 65 mm Locking Hexalobe Screw	30-0253