

Item Catalogue

NUPCO CODE: 4232150013000 | SN: 690

NUPCO ITEM DESCRIPTION

SCREW, CORTICAL 3.5MM OR EQUIVALENT FOR TITANIUM HOOK CLAVICLE PLATE, ASSORTED SIZES

VENDOR PRODUCT DESCRIPTION

Acumed Clavicle Hook Plate 3.5 mm cortical/nonlocking screw families

SKU CODES INCLUDED FOR THIS NUPCO CODE

- 30-0255
- 30-0256
- 30-0257
- 30-0258
- 30-0259
- 30-0260
- 30-0261
- 30-0262
- 30-0263
- 30-0264
- 30-0265
- 30-0266
- 30-0267
- 30-0268
- 30-0269
- 30-0270
- 30-0271
- 30-3542
- 30-3544
- 30-3546
- 30-3548
- 30-0273
- 30-0274
- 30-0275
- 30-0276

Surgical Technique



Ordering Information [continued]

Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264

3.5 mm (Nonlocking) Cortical (Hex) Screws

3.5 mm x 8 mm Cortical Screw	CO-3080
3.5 mm x 10 mm Cortical Screw	CO-3100
3.5 mm x 12 mm Cortical Screw	CO-3120
3.5 mm x 14 mm Cortical Screw	CO-3140
3.5 mm x 16 mm Cortical Screw	CO-3160
3.5 mm x 18 mm Cortical Screw	CO-3180
3.5 mm x 20 mm Cortical Screw	CO-3200
3.5 mm x 22 mm Cortical Screw	CO-3220
3.5 mm x 24 mm Cortical Screw	CO-3240
3.5 mm x 26 mm Cortical Screw	CO-3260

3.5 mm Locking Cortical (Hex) Screws

3.5 mm x 8 mm Locking Cortical Screw	COL-3080
3.5 mm x 10 mm Locking Cortical Screw	COL-3100
3.5 mm x 12 mm Locking Cortical Screw	COL-3120
3.5 mm x 14 mm Locking Cortical Screw	COL-3140
3.5 mm x 16 mm Locking Cortical Screw	COL-3160
3.5 mm x 18 mm Locking Cortical Screw	COL-3180
3.5 mm x 20 mm Locking Cortical Screw	COL-3200
3.5 mm x 22 mm Locking Cortical Screw	COL-3220
3.5 mm x 24 mm Locking Cortical Screw	COL-3240
3.5 mm x 26 mm Locking Cortical Screw	COL-3260

Ordering Information [continued]

Sterile Screws

Sterile 3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232-S
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233-S
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234-S
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235-S
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236-S
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237-S
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238-S
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239-S
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240-S
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241-S

Sterile 3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255-S
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256-S
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257-S
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258-S
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259-S
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260-S
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261-S
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262-S
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263-S
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264-S

Sterile 3.5 mm (Nonlocking) Cortical (Hex) Screws


3.5 mm x 8 mm Cortical Screw	CO-3080-S
3.5 mm x 10 mm Cortical Screw	CO-3100-S
3.5 mm x 12 mm Cortical Screw	CO-3120-S
3.5 mm x 14 mm Cortical Screw	CO-3140-S
3.5 mm x 16 mm Cortical Screw	CO-3160-S
3.5 mm x 18 mm Cortical Screw	CO-3180-S
3.5 mm x 20 mm Cortical Screw	CO-3200-S
3.5 mm x 22 mm Cortical Screw	CO-3220-S
3.5 mm x 24 mm Cortical Screw	CO-3240-S
3.5 mm x 26 mm Cortical Screw	CO-3260-S

Sterile 3.5 mm Locking Cortical (Hex) Screws

3.5 mm x 8 mm Locking Cortical Screw	COL-3080-S
3.5 mm x 10 mm Locking Cortical Screw	COL-3100-S
3.5 mm x 12 mm Locking Cortical Screw	COL-3120-S
3.5 mm x 14 mm Locking Cortical Screw	COL-3140-S
3.5 mm x 16 mm Locking Cortical Screw	COL-3160-S
3.5 mm x 18 mm Locking Cortical Screw	COL-3180-S
3.5 mm x 20 mm Locking Cortical Screw	COL-3200-S
3.5 mm x 22 mm Locking Cortical Screw	COL-3220-S
3.5 mm x 24 mm Locking Cortical Screw	COL-3240-S
3.5 mm x 26 mm Locking Cortical Screw	COL-3260-S

Note: To learn more about the full line of Acumed innovative surgical solutions, please contact your authorized Acumed distributor, call 888.627.9957, or visit www.acumed.net.

Surgical Technique



Superior Distal Clavicle Plate Surgical Technique [continued]

4 Plate Placement

Once the plate's ideal positioning has been selected, it is provisionally stabilized to the clavicle with .045" x 6" ST Guide Wire (WS-1106ST) or .059" x 5" ST Guide Wire (WS-1505ST) K-wires. Under radiographic evaluation, the most lateral K-wire hole of each Superior Distal Clavicle Plate (70-0XXX) affords the opportunity to verify that the placement of the screws will not protrude into the AC joint by inserting a K-wire to confirm plate placement.

Note: The Reduction Forceps With Serrated Jaw (PL-CL04) should be used for plate placement and are not designed to be used to reduce the plate to the bone or to hold the plate while attempting to bend or contour it to match the patient's anatomy. Plate Benders (PL-204X) are available in the event that plate contouring is required to achieve the desired fit to the clavicle.


Caution: If bending of the plate is necessary, please observe the following:

- ▶ Do not bend plates more than 30 degrees
- ▶ Bend radii should be greater than 1 inch
- ▶ Do not bend more than once
- ▶ Avoid bending across locking holes

Figure 5



Superior Distal Clavicle Plate Surgical Technique [continued]



5 Nonlocking Screw Insertion

For early stability, place the first two screws medial and lateral to the fracture site. If bicortical screws are used, take precautions to avoid over-penetration of the inferior cortex. Place the Clavicle Retractor (PL-CL03) under the inferior surface of the clavicle to protect the neurovascular structures from overpenetration when drilling.

Caution: Replace the drill if it comes in contact with the clavicle retractor.

If drilling for 3.5 mm size screws, the Offset Drill Guide (PL-2095) is available to assist in providing compression. The green side of the drill guide is neutral, while the gold side will generate compression when used in a compression slot on the plate.

Screw Size Option	Drill	Depth Gauge	Driver	Driver Handle
3.5 mm Nonlocking Hexalobe Screw (30-02XX)	2.8 mm Quick Release Drill (80-0387)	Depth Gauge 6–65 mm (80-0623)	T15 Stick Fit Hexalobe Driver (80-0760)	Large Cannulated Quick Release Driver Handle (MS-3200)
3.0 mm Nonlocking Hexalobe Screw (30-03XX)	2.3 mm Quick Release Drill (80-0627)	Depth Gauge 6–65 mm (80-0623)	T15 Stick Fit Hexalobe Driver (80-0760)	Large Cannulated Quick Release Driver Handle (MS-3200)
Optional 3.5 mm (Nonlocking) Cortical (Hex) Screw (CO-3XX0)	2.8 mm x 5" Quick Release Drill (MS-DC28)	6–70 mm Depth Gauge, 2 mm Increments (MS-9022)	2.5 mm Quick Release Hex Driver (HPC-0025)	Large Cannulated Quick Release Driver Handle (MS-3200)
Optional 2.7 mm (Nonlocking) Cortical (Hex) Screw (CO-27XX)	2.0 mm x 5" Quick Release Drill (MS-DC5020)	6–70 mm Depth Gauge, 2 mm Increments (MS-9022)	2.5 mm Quick Release Hex Driver (HPC-0025)	Large Cannulated Quick Release Driver Handle (MS-3200)

Ordering Information [continued]

Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264

Hexalobe Short Screw Caddies

3.5 mm Hexalobe Short Screw Caddy	80-0843
3.5 mm Hexalobe Short Screw Caddy Lid	80-0856
3.0 mm Hexalobe Short Screw Caddy	80-1066
3.0 mm Hexalobe Short Screw Caddy Lid	80-1067

3.0 mm Locking Hexalobe Screws

3.0 mm x 8 mm Locking Hexalobe Screw	30-0278
3.0 mm x 10 mm Locking Hexalobe Screw	30-0279
3.0 mm x 12 mm Locking Hexalobe Screw	30-0280
3.0 mm x 14 mm Locking Hexalobe Screw	30-0281
3.0 mm x 16 mm Locking Hexalobe Screw	30-0282
3.0 mm x 18 mm Locking Hexalobe Screw	30-0283
3.0 mm x 20 mm Locking Hexalobe Screw	30-0284
3.0 mm x 22 mm Locking Hexalobe Screw	30-0285
3.0 mm x 24 mm Locking Hexalobe Screw	30-0286
3.0 mm x 26 mm Locking Hexalobe Screw	30-0287

3.0 mm Nonlocking Hexalobe Screws

3.0 mm x 8 mm Nonlocking Hexalobe Screw	30-0301
3.0 mm x 10 mm Nonlocking Hexalobe Screw	30-0302
3.0 mm x 12 mm Nonlocking Hexalobe Screw	30-0303
3.0 mm x 14 mm Nonlocking Hexalobe Screw	30-0304
3.0 mm x 16 mm Nonlocking Hexalobe Screw	30-0305
3.0 mm x 18 mm Nonlocking Hexalobe Screw	30-0306
3.0 mm x 20 mm Nonlocking Hexalobe Screw	30-0307
3.0 mm x 22 mm Nonlocking Hexalobe Screw	30-0308
3.0 mm x 24 mm Nonlocking Hexalobe Screw	30-0309
3.0 mm x 26 mm Nonlocking Hexalobe Screw	30-0310

Ordering Information [continued]

Sterile Screws

3.5 mm Locking Hexalobe Screws

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232-S
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233-S
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234-S
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235-S
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236-S
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237-S
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238-S
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239-S
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240-S
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241-S

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255-S
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256-S
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257-S
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258-S
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259-S
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260-S
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261-S
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262-S
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263-S
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264-S

3.0 mm Locking Hexalobe Screws

3.0 mm x 8 mm Locking Hexalobe Screw	30-0278-S
3.0 mm x 10 mm Locking Hexalobe Screw	30-0279-S
3.0 mm x 12 mm Locking Hexalobe Screw	30-0280-S
3.0 mm x 14 mm Locking Hexalobe Screw	30-0281-S
3.0 mm x 16 mm Locking Hexalobe Screw	30-0282-S
3.0 mm x 18 mm Locking Hexalobe Screw	30-0283-S
3.0 mm x 20 mm Locking Hexalobe Screw	30-0284-S
3.0 mm x 22 mm Locking Hexalobe Screw	30-0285-S
3.0 mm x 24 mm Locking Hexalobe Screw	30-0286-S
3.0 mm x 26 mm Locking Hexalobe Screw	30-0287-S

3.0 mm Nonlocking Hexalobe Screws

3.0 mm x 8 mm Nonlocking Hexalobe Screw	30-0301-S
3.0 mm x 10 mm Nonlocking Hexalobe Screw	30-0302-S
3.0 mm x 12 mm Nonlocking Hexalobe Screw	30-0303-S
3.0 mm x 14 mm Nonlocking Hexalobe Screw	30-0304-S
3.0 mm x 16 mm Nonlocking Hexalobe Screw	30-0305-S
3.0 mm x 18 mm Nonlocking Hexalobe Screw	30-0306-S
3.0 mm x 20 mm Nonlocking Hexalobe Screw	30-0307-S
3.0 mm x 22 mm Nonlocking Hexalobe Screw	30-0308-S
3.0 mm x 24 mm Nonlocking Hexalobe Screw	30-0309-S
3.0 mm x 26 mm Nonlocking Hexalobe Screw	30-0310-S

 acumed[®]
Small Fragment Base Set



Acumed System
Compatibility Indicator



Acumed Variable
Angle Screw Compatible

Surgical Technique



A COLSON ASSOCIATE

Variable Angle Screw Surgical Technique [continued]

4 Insert Variable Angle Screw

Note: Final tightening of the 2.7 mm and 3.5 mm variable angle hexalobe screws must be done manually and not under power. The Torque Limiting Quick Connect ensures a consistent insertion torque to provide a uniform screw-plate interface and may prevent overtightening of the screw. See below for the torque limit for each screw.

To insert a 2.7 mm variable angle hexalobe screw, assemble the 1.70 N•m Torque Limiting Quick Connect (80-2366) to the Handle for Torque Limiting Quick Connect (80-2368). Connect the T8 Stick Fit Hexalobe Driver (80-0759) to the Torque Limiting Quick Connect assembly.

To insert a 3.5 mm variable angle hexalobe screw, assemble the 2.26 N•m Torque Limiting Quick Connect (80-2367) to the Handle for Torque Limiting Quick Connect (80-2368). Connect the T15 Stick Fit Hexalobe Driver (80-0760) to the Torque Limiting Quick Connect assembly.

Advance the screw by hand until achieving an audible click and/or tactile feedback. Final tightening should be completed with the Torque Limiting Handle, which is designed to provide a secure lock between the plate and screw. Upon final seating, confirm proper screw placement and screw length under fluoroscopy.

Note: Do not use a Torque Limiting Quick Connect for screw removal.

Caution: Use the maximum number of screws based on the indication to reduce the risk of screw breakage during healing.




Figure 18

Screw	Torque Limit	Color Band
2.7 mm Variable Angle Hexalobe Screw (3013-27XXX)	1.70 N•m	Brown
3.5 mm Variable Angle Hexalobe Screw (3013-35XXX)	2.26 N•m	Black



Handle for Torque Limiting Quick Connect (80-2368)



T8 Stick Fit Hexalobe Driver (80-0759)



T15 Stick Fit Hexalobe Driver (80-0760)



2.7 mm Variable Angle Hexalobe Screw, 10–60 mm lengths (3013-27XXX)



3.5 mm Variable Angle Hexalobe Screw, 10–65 mm lengths (3013-35XXX)



1.70 N•m Torque Limiting Quick Connect (80-2366)



2.26 N•m Torque Limiting Quick Connect (80-2367)

Ordering Information

Tray Components

Implants*

1	Fragment Plate 2.7 mm	7010-0106N	10	One-Third Tubular Plate 3-Hole	7008-0103
2	L Fragment Plate 2.7 mm Left	7010-0107L	11	One-Third Tubular Plate 4-Hole	7008-0104
3	L Fragment Plate 2.7 mm Right	7010-0107R	12	One-Third Tubular Plate 5-Hole	7008-0105
4	T Fragment Plate 2.7 mm	7010-0108N	13	One-Third Tubular Plate 6-Hole	7008-0106
5	20 mm AcuTwist® Acutrak Compression Screw	AI-0020	14	One-Third Tubular Plate 7-Hole	7008-0107
6	26 mm AcuTwist® Acutrak Compression Screw	AI-0026	15	One-Third Tubular Plate 8-Hole	7008-0108
7	30 mm AcuTwist® Acutrak Compression Screw	AI-0030	16	One-Third Tubular Plate 10-Hole	7008-0110
8	70.0 mm Tension Band Pin (1.6 mm diameter)	30-0098	17	One-Third Tubular Plate 12-Hole	7008-0112
9	90.0 mm Tension Band Pin	30-0099	18	Cannulated Screw Washer 7.0 mm OD x 3.6 mm ID	7003-07036

Instruments

19	30 mm AcuTwist® Acutrak® Tap	AI-NG30	28	Periosteal Elevator	MS-46212
20	Plate Bender, Large	PL-2045	29	8 mm Hohman Retractor	PL-CL05
21	Tension Band Pin Snapper	80-0411	30	15 mm Hohman Retractor	MS-46827
22	Needle Nose Pliers, 5.5	MS-48245	31	Sharp Hook	PL-CL06
23	AcuTwist® Acutrak® Compression Screw Extractor	AI-EX20	32	Pointed Forceps w/ Ratchet, Wide Long	80-2375
24	Fragment Plate Cutter	80-2380	33	Pointed Forceps w/ Ratchet, Narrow Long	80-2376
25	Fragment Plate Bender, Long	80-2381	34	Reduction Forceps w/ Ratchet, Long	80-2377
26	Fragment Plate Bender, Short	80-2382	35	Large Screw Holding Forceps	MS-45210
27	Polarus® 3 Reduction Device	80-1601			

*Implants and screws are also available sterile-packed. Add an “-S” at end of product number for sterile product. For more details on sterile products, including pricing, contact our Business Services Department toll free at 888.627.9957.

Ordering Information [continued]

Screws

4 3.5 mm Nonlocking Hexalobe Screws*

3.5 mm x 8 mm Nonlocking Hexalobe Screw	30-0255
3.5 mm x 10 mm Nonlocking Hexalobe Screw	30-0256
3.5 mm x 12 mm Nonlocking Hexalobe Screw	30-0257
3.5 mm x 14 mm Nonlocking Hexalobe Screw	30-0258
3.5 mm x 16 mm Nonlocking Hexalobe Screw	30-0259
3.5 mm x 18 mm Nonlocking Hexalobe Screw	30-0260
3.5 mm x 20 mm Nonlocking Hexalobe Screw	30-0261
3.5 mm x 22 mm Nonlocking Hexalobe Screw	30-0262
3.5 mm x 24 mm Nonlocking Hexalobe Screw	30-0263
3.5 mm x 26 mm Nonlocking Hexalobe Screw	30-0264
3.5 mm x 28 mm Nonlocking Hexalobe Screw	30-0265
3.5 mm x 30 mm Nonlocking Hexalobe Screw	30-0266
3.5 mm x 32 mm Nonlocking Hexalobe Screw	30-0267
3.5 mm x 34 mm Nonlocking Hexalobe Screw	30-0268
3.5 mm x 36 mm Nonlocking Hexalobe Screw	30-0269
3.5 mm x 38 mm Nonlocking Hexalobe Screw	30-0270
3.5 mm x 40 mm Nonlocking Hexalobe Screw	30-0271
3.5 mm x 42 mm Nonlocking Hexalobe Screw	30-3542
3.5 mm x 44 mm Nonlocking Hexalobe Screw	30-3544
3.5 mm x 46 mm Nonlocking Hexalobe Screw	30-3546
3.5 mm x 48 mm Nonlocking Hexalobe Screw	30-3548
3.5 mm x 50 mm Nonlocking Hexalobe Screw	30-0273
3.5 mm x 55 mm Nonlocking Hexalobe Screw	30-0274
3.5 mm x 60 mm Nonlocking Hexalobe Screw	30-0275
3.5 mm x 65 mm Nonlocking Hexalobe Screw	30-0276

5 3.5 mm Locking Hexalobe Screws*

3.5 mm x 8 mm Locking Hexalobe Screw	30-0232
3.5 mm x 10 mm Locking Hexalobe Screw	30-0233
3.5 mm x 12 mm Locking Hexalobe Screw	30-0234
3.5 mm x 14 mm Locking Hexalobe Screw	30-0235
3.5 mm x 16 mm Locking Hexalobe Screw	30-0236
3.5 mm x 18 mm Locking Hexalobe Screw	30-0237
3.5 mm x 20 mm Locking Hexalobe Screw	30-0238
3.5 mm x 22 mm Locking Hexalobe Screw	30-0239
3.5 mm x 24 mm Locking Hexalobe Screw	30-0240
3.5 mm x 26 mm Locking Hexalobe Screw	30-0241
3.5 mm x 28 mm Locking Hexalobe Screw	30-0242
3.5 mm x 30 mm Locking Hexalobe Screw	30-0243
3.5 mm x 32 mm Locking Hexalobe Screw	30-0244
3.5 mm x 34 mm Locking Hexalobe Screw	30-0245
3.5 mm x 36 mm Locking Hexalobe Screw	30-0246
3.5 mm x 38 mm Locking Hexalobe Screw	30-0247
3.5 mm x 40 mm Locking Hexalobe Screw	30-0248
3.5 mm x 42 mm Locking Hexalobe Screw	30-3642
3.5 mm x 44 mm Locking Hexalobe Screw	30-3644
3.5 mm x 46 mm Locking Hexalobe Screw	30-3646
3.5 mm x 48 mm Locking Hexalobe Screw	30-3648
3.5 mm x 50 mm Locking Hexalobe Screw	30-0250
3.5 mm x 55 mm Locking Hexalobe Screw	30-0251
3.5 mm x 60 mm Locking Hexalobe Screw	30-0252
3.5 mm x 65 mm Locking Hexalobe Screw	30-0253