

Item Catalogue

NUPCO CODE: 4232150015900 | SN: 799

NUPCO ITEM DESCRIPTION

PLATE, LCP TITANIUM FOR FOREFOOT AND MIDFOOT, OR EQUIVALENTS, ASSORTED SIZES


VENDOR PRODUCT DESCRIPTION

Forefoot/Midfoot plate families listed in source catalogue

SKU CODES INCLUDED FOR THIS NUPCO CODE

- 70-0141
- 70-0142
- 70-0012
- 70-0013
- 70-0143
- 70-0144
- 70-0109
- 70-0110
- 70-0329
- 70-0330
- 70-0018
- 70-0019
- 70-0036
- 70-0037
- 70-0010
- 70-0011
- 70-0323
- 70-0324
- 70-0007
- 70-0008
- 70-0009
- 70-0049
- 70-0050
- 70-0051

Surgical Technique



System Features

Plate Family

Multiple Plate Options offer surgeons choices in forefoot/midfoot plating from an expanding line of MTP, TMT, and Osteotomy Plates.

Precontoured Plates are designed to match the anatomy of the patient. Based on multiple cadaveric and clinical trials, the MTP Plates are precontoured with both 4 degrees and 9 degrees of dorsiflexion and 10 degrees of lateral translation to help restore the functional angle of the MTP joint after fusion.


Blue plates are left specific, green plates are right specific, turquoise plates are neither right nor left specific.



Plate options include:

Metatarsophalangeal (MTP) Plates (available in both 4° and 9° dorsiflexion)

- ▶ Standard MTP Plate
- ▶ Dorsal MTP Plate
- ▶ MTP Revision Plate
- ▶ Petite MTP Plate (Available in 4° dorsiflexion only)




System Features [continued]

Forefoot and Midfoot Plates



Osteotomy Plates

- ▶ Osteotomy Plate
- ▶ Osteotomy Plate With Compression Slot




Locking Proximal
MT Wedge, Left
(70-00XX)

Locking Proximal
MT Wedge, Right
(70-00XX)

Tarsometatarsal (TMT) Plates

- ▶ First Ray TMT Plate
- ▶ Second/Third Ray TMT Plate
- ▶ Lapidus Plate
- ▶ Extended First Ray TMT Plate
- ▶ Extended Second/Third Ray TMT Plate
- ▶ Extended Lapidus Plate



Ordering Information

Tray Components

Forefoot/Midfoot Plates

1	Locking Dorsal MTP Plate, 4 Degrees, Left	70-0141	17	Locking Proximal MT Wedge Plate, 0 Degrees, Left	70-0323
2	Locking Dorsal MTP Plate, 4 Degrees, Right	70-0142	18	Locking Proximal MT Wedge Plate, 0 Degrees, Right	70-0324
3	Locking Dorsal MTP Fusion Plate, Left	70-0012	19	5-Hole Locking 1st TMT Plate	70-0007
4	Locking Dorsal MTP Fusion Plate, Right	70-0013	20	4-Hole Locking 1st TMT Plate	70-0008
5	Locking MTP Revision Plate, 4 Degrees, Left	70-0143	21	4-Hole Locking 2nd & 3rd TMT Plate	70-0009
6	Locking MTP Revision Plate, 4 Degrees, Right	70-0144	22	8-Hole Locking 1st TMT Plate	70-0049
7	Locking MTP/MPJ Combo Plate 4 Degrees, Left	70-0109	23	7-Hole Locking 1st TMT Plate	70-0050
8	Locking MTP/MPJ Combo Plate 4 Degrees, Right	70-0110	24	7-Hole Locking 2nd & 3rd TMT Plate	70-0051
9	Locking MPJ Petite Plate 4 Degrees, Left	70-0329	25	MTP Reamer Concave, 16 mm with QR	80-0569
10	Locking MPJ Petite Plate 4 Degrees, Right	70-0330	26	MTP Reamer Convex, 16 mm with QR	80-0570
11	Locking MTP Revision Fusion Plate, Left	70-0018	27	MTP Reamer Concave, 18 mm with QR	80-0571
12	Locking MTP Revision Fusion Plate, Right	70-0019	28	MTP Reamer Convex, 18 mm with QR	80-0572
13	Locking MTP/MPJ Combo Fusion Plate, Left	70-0036	29	MTP Reamer Concave, 20 mm with QR	80-0573
14	Locking MTP/MPJ Combo Fusion Plate, Right	70-0037	30	MTP Reamer Convex, 20 mm with QR	80-0574
15	Locking Proximal MT Wedge Plate, Left	70-0010	31	MTP Reamer Concave, 22 mm with QR	80-0575
16	Locking Proximal MT Wedge Plate, Right	70-0011	32	MTP Reamer Convex, 22 mm with QR	80-0576

