

Item Catalogue

NUPCO CODE: 4232170000100 | SN: 835

NUPCO ITEM DESCRIPTION

SCREW, CORTICAL FOR LOW PROFILE PELVIS RECONSTRUCTION PLATING, ASSORTED SIZES

VENDOR PRODUCT DESCRIPTION


Pelvic nonlocking/column screw options listed in source catalogue

SKU CODES INCLUDED FOR THIS NUPCO CODE

30-0344	30-0881
30-0345	30-0882
30-0346	30-0883
30-0347	30-0884
30-0348	30-0885
30-0349	30-0886
30-0350	30-0887
30-0351	30-0888
30-0352	30-0889
30-0353	30-0890
30-0354	30-0891
30-0355	30-0892
30-0256	30-0893
30-0257	7003-10046
30-0258	30-0894
30-0259	30-0895
30-0260	30-0896
30-0261	30-0897
30-0262	30-0898
30-0263	30-0899
30-0264	30-0900
30-0265	30-0901
30-0266	30-0902
30-0267	30-0903
30-0268	30-0904
30-0269	30-0905
30-0270	30-0906
30-0271	30-0907
30-0272	30-0908
30-0273	30-0909
30-0274	30-0910
30-0275	30-0911
30-0877	30-0912
30-0878	30-0913
30-0879	30-0914
30-0880	



Surgical Technique



4.3 mm Hexalobe Column Screw Surgical Technique [continued]

3 Screw Insertion

Connect the T15 Stick Fit Hexalobe Driver, 6 inch (80-2087) to the Small Ratchet Handle with Quick Release Connection (80-0398) and insert a 4.3 mm Hexalobe Column Screw (30-0XXX). Fully seat the screw.

Note: In cases of dense bone, it may be necessary to open the near cortex using the 4.3 mm Quick Release Drill, Long (80-1133).




Figure 77

4 Screw Insertion

Repeat Steps 1–3 at the surgeon's discretion. Continue the definitive fixation of the fracture using the appropriate surgical technique previously described. Confirm all screws have been seated. Close the surgical site(s) using preferred techniques.

Note: Optional Cannulated Screw Washers (7003-10046) are included in the Acumed Cannulated Screw System, and may be used with the 4.3 mm Hexalobe Column Screws (30-0XXX) at the surgeon's discretion.




Figure 78

5 Optional: Removal Instructions

4.3 mm Hexalobe Column Screws may be removed using either the T15 Hexalobe Driver Tip, Long (80-1561), T15 Stick Fit Hexalobe Driver, 6 Inch (80-2087), or T15 Stick Fit Hexalobe Driver, 9 Inch (80-2088) and the Small Ratchet Handle with Quick Release Connection (80-0398).

Ordering Information [continued]

Screws

2.7 mm Nonlocking Hexalobe Screws

2.7 mm x 10 mm Nonlocking Hexalobe Screw*	30-0344	2.7 mm x 22 mm Nonlocking Hexalobe Screw*	30-0350
2.7 mm x 12 mm Nonlocking Hexalobe Screw*	30-0345	2.7 mm x 24 mm Nonlocking Hexalobe Screw*	30-0351
2.7 mm x 14 mm Nonlocking Hexalobe Screw*	30-0346	2.7 mm x 26 mm Nonlocking Hexalobe Screw*	30-0352
2.7 mm x 16 mm Nonlocking Hexalobe Screw*	30-0347	2.7 mm x 28 mm Nonlocking Hexalobe Screw*	30-0353
2.7 mm x 18 mm Nonlocking Hexalobe Screw*	30-0348	2.7 mm x 30 mm Nonlocking Hexalobe Screw*	30-0354
2.7 mm x 20 mm Nonlocking Hexalobe Screw*	30-0349	2.7 mm x 32 mm Nonlocking Hexalobe Screw*	30-0355

3.5 mm Nonlocking Hexalobe Screws

3.5 mm x 10 mm Nonlocking Hexalobe Screw*	30-0256	3.5 mm x 60 mm Nonlocking Hexalobe Screw*	30-0275
3.5 mm x 12 mm Nonlocking Hexalobe Screw*	30-0257	3.5 mm x 65 mm Nonlocking Hexalobe Screw*	30-0276
3.5 mm x 14 mm Nonlocking Hexalobe Screw*	30-0258	3.5 mm x 70 mm Nonlocking Hexalobe Screw*	30-0877
3.5 mm x 16 mm Nonlocking Hexalobe Screw*	30-0259	3.5 mm x 75 mm Nonlocking Hexalobe Screw*	30-0878
3.5 mm x 18 mm Nonlocking Hexalobe Screw*	30-0260	3.5 mm x 80 mm Nonlocking Hexalobe Screw*	30-0879
3.5 mm x 20 mm Nonlocking Hexalobe Screw*	30-0261	3.5 mm x 85 mm Nonlocking Hexalobe Screw*	30-0880
3.5 mm x 22 mm Nonlocking Hexalobe Screw*	30-0262	3.5 mm x 90 mm Nonlocking Hexalobe Screw*	30-0881
3.5 mm x 24 mm Nonlocking Hexalobe Screw*	30-0263	3.5 mm x 95 mm Nonlocking Hexalobe Screw*	30-0882
3.5 mm x 26 mm Nonlocking Hexalobe Screw*	30-0264	3.5 mm x 100 mm Nonlocking Hexalobe Screw*	30-0883
3.5 mm x 28 mm Nonlocking Hexalobe Screw*	30-0265	3.5 mm x 105 mm Nonlocking Hexalobe Screw*	30-0884
3.5 mm x 30 mm Nonlocking Hexalobe Screw*	30-0266	3.5 mm x 110 mm Nonlocking Hexalobe Screw*	30-0885
3.5 mm x 32 mm Nonlocking Hexalobe Screw*	30-0267	3.5 mm x 115 mm Nonlocking Hexalobe Screw*	30-0886
3.5 mm x 34 mm Nonlocking Hexalobe Screw*	30-0268	3.5 mm x 120 mm Nonlocking Hexalobe Screw*	30-0887
3.5 mm x 36 mm Nonlocking Hexalobe Screw*	30-0269	3.5 mm x 125 mm Nonlocking Hexalobe Screw*	30-0888
3.5 mm x 38 mm Nonlocking Hexalobe Screw*	30-0270	3.5 mm x 130 mm Nonlocking Hexalobe Screw*	30-0889
3.5 mm x 40 mm Nonlocking Hexalobe Screw*	30-0271	3.5 mm x 135 mm Nonlocking Hexalobe Screw*	30-0890
3.5 mm x 45 mm Nonlocking Hexalobe Screw*	30-0272	3.5 mm x 140 mm Nonlocking Hexalobe Screw*	30-0891
3.5 mm x 50 mm Nonlocking Hexalobe Screw*	30-0273	3.5 mm x 145 mm Nonlocking Hexalobe Screw*	30-0892
3.5 mm x 55 mm Nonlocking Hexalobe Screw*	30-0274	3.5 mm x 150 mm Nonlocking Hexalobe Screw*	30-0893

*Implants and screws are also available sterile-packed. Add an "-S" at end of product number for sterile product. For more details on sterile products, including pricing, contact our Business Services Department toll free at 888.627.9957.

Ordering Information [continued]

Screws	Optional
4.3 mm Hexalobe Column Screws	Cannulated Screw Washers 10.0 mm OD x 4.6 mm ID* 7003-10046
4.3 mm x 50 mm Hexalobe Column Screw*	30-0894
4.3 mm x 55 mm Hexalobe Column Screw*	30-0895
4.3 mm x 60 mm Hexalobe Column Screw*	30-0896
4.3 mm x 65 mm Hexalobe Column Screw*	30-0897
4.3 mm x 70 mm Hexalobe Column Screw*	30-0898
4.3 mm x 75 mm Hexalobe Column Screw*	30-0899
4.3 mm x 80 mm Hexalobe Column Screw*	30-0900
4.3 mm x 85 mm Hexalobe Column Screw*	30-0901
4.3 mm x 90 mm Hexalobe Column Screw*	30-0902
4.3 mm x 95 mm Hexalobe Column Screw*	30-0903
4.3 mm x 100 mm Hexalobe Column Screw*	30-0904
4.3 mm x 105 mm Hexalobe Column Screw*	30-0905
4.3 mm x 110 mm Hexalobe Column Screw*	30-0906
4.3 mm x 115 mm Hexalobe Column Screw*	30-0907
4.3 mm x 120 mm Hexalobe Column Screw*	30-0908
4.3 mm x 125 mm Hexalobe Column Screw*	30-0909
4.3 mm x 130 mm Hexalobe Column Screw*	30-0910
4.3 mm x 135 mm Hexalobe Column Screw*	30-0911
4.3 mm x 140 mm Hexalobe Column Screw*	30-0912
4.3 mm x 145 mm Hexalobe Column Screw*	30-0913
4.3 mm x 150 mm Hexalobe Column Screw*	30-0914

*Implants and screws are also available sterile-packed. Add an “-S” at end of product number for sterile product.
For more details on sterile products, including pricing, contact our Business Services Department toll free at 888.627.9957.